

# Q-sep™ QuEChERS Products

## BUY 1, GET 1 FREE!

extraction salt packet products



### Introducing New Q-sep Extraction Salt Packets!

- Convenient and easy to use.
- Eliminates the need for a second tube to transfer salts.

*Supports AOAC (2007.01) & European (EN15662) QuEChERS Methods*

**Buy 1, get 1 FREE**  
offer applies to  
the following  
products:

Description	Material	Methods	qty.	cat#
Q-sep Q110 kit	4g MgSO <sub>4</sub> , 1g NaCl, 1g TSCD, 0.5g DHS with 50mL Centrifuge Tube	European EN 15662	50 packets & 50 tubes	26235-408
Q-sep Q110 packets	4g MgSO <sub>4</sub> , 1g NaCl, 1g TSCD, 0.5g DHS	European EN 15662	50 packets	26236-408
Q-sep Q150 kit	6g MgSO <sub>4</sub> , 1.5g NaOAc, 1g TSCD, 0.5g DHS with 50mL Centrifuge Tube	AOAC 2007.01	50 packets & 50 tubes	26237-408
Q-sep Q150 packets	6g MgSO <sub>4</sub> , 1.5g NaOAc	AOAC 2007.01	50 packets	26238-408

Offer valid July 1—October 31, 2010.

TSCD = trisodium citrate dihydrate  
DHS = disodium hydrogen citrate sesquihydrate  
NaOAc = sodium acetate

Lit. Cat.# FFTS1281  
© 2010 Restek Corporation.

ISO 9001:2008  
cert.# FM80397

# RESTEK

## Chromatography Products

www.restek.com